

ComfortLine with rechargeable battery: ...more independence

Besides the well known advantages of the ComfortLine, for example the **5 easy to access motors**, the **multifunctional armrests** and the activation of the **shock position by foot pedal** the ComfortLine with rechargeable battery provides independence from the electric supply. Thus the therapy chair is ready for operation and can be moved freely at any time.

Different colour applications make the chair fit perfectly into the desired environment.

ComfortLine with rechargeable battery and bi-colour design



Standard colour:

Chrom grey



Special colours for
bi-colour applications:

Feuer



Atoll



Lima



Amazonas



Technical Data

Dimensions

Seat section width	60 cm
Chair width	90 cm
Total length	215 cm
Total weight	98 kg

Adjustment features

Back section	75° - 0°
Seat section	30° - 5°
Leg section	5° - 30°
Foot rest electrical	25 cm
Foot rest mechanical additional	18 cm
Height adjustment	60 - 80 cm

Safety

Safe maximum load	180 kg
Angle of shock position	- 12°
Time for shock position	< 20 sec
Protection class	2
Protection category	IPX 4
In compliance with	EN 60601
Motors	5

 **ComfortLine®**

Therapy chair
with rechargeable battery

Dialysis
Transfusion
Oncology




bionic

Bionic Medizintechnik GmbH ... a JMS-corporation · Max-Planck-Straße 21
D-61381 Friedrichsdorf/Germany · Phone +49 (0)6172-7576-0
Fax +49 (0)6172-7576-10 · info@bionic-jms.com · www.bionic-jms.com

Right for technical modifications reserved
B-MORAKKE1


bionic



Therapy chair with rechargeable battery



Characteristics



The rechargeable battery guarantees independence from the electric supply and gives the possibility to move the therapy chair freely. The optional second rechargeable battery provides greater independence.

Rechargeable battery - dimensions

Width	9.9 cm
Length	47.3 cm
Height	12.1 cm
Weight	6.2 kg

Electrical specifications

Battery charger

Voltage	230 V 50/60 Hz
Current	830 mA
Primary power	26 W
Secondary power	20 W

Rechargeable battery

Voltage	24 V
Power	7.2 Ah
Charge time	10 - 12 h